

genesig[®] SARS-CoV-2 Winterplex



Respiratory viruses, including COVID-19, are known to have a high incidence of infection during winter, especially in temperate regions. However, these winter respiratory illnesses have many similar signs and symptoms. So how do we know if it's COVID-19, Flu A, Flu B or RSV? We have just the right multiplex assay for you.

Introducing the next generation genesig[®] SARS-CoV-2 Winterplex, a multiplex assay that enables simultaneous real-time PCR detection of multiple winter viruses in one kit, successfully distinguishing SARS CoV-2 and all 20/21 vaccine strains of Flu A and Flu B, and RSV in as quickly as 90 minutes.

DETECTION PROFILE

- SARS-CoV-2 (ORF1ab, S and M gene targets)
- Influenza A
- Influenza B
- Respiratory Syncytial Virus (RSV A and B)

KEY FEATURES

- Precise reproducible results
- Distinguish SARS-CoV-2, Flu A, Flu B and RSV (including all the latest strains) with one assay in as quickly as 90 minutes
- Kit includes Lyophilised qPCR Master Mix for ambient shipping
- High sensitivity and accuracy

ORDERING INFORMATION

Product code: D00020
Pack size: 96 reactions (Master Mix included)

PLEASE CONTACT ON: +44 (0) 2380 748 830

Novacyt Group (UK)
Unit 1, Watchmoor Point, Watchmoor Road
Camberley, Surrey GU15 3AD, United Kingdom

www.novacyt.com



NOVACYT
GROUP