



## VersaLab™ Mobile Processing Laboratories

**VersaLab™** aspires to support private sector testing of infectious diseases, initially focused on COVID-19. The VersaLab™ service includes the VersaLab™ mobile processing laboratories (MPL) to provide rapid turnaround polymerase chain reaction (PCR) testing on-site. VersaLab™ MPLs offer the private sectors the opportunities to take control of COVID-19 testing to keep businesses open and running efficiently.

### Key Benefits

- PROmate™ and genesig® PCR instruments work seamlessly in a streamlined workflow making lab testing more efficient and cost effective.
- Results within 80 minutes, without the need of a CL2 Lab.
- Quick and easy sample processing without requirement for a highly trained technician
- Accommodating at least two workstations
- Can process up to 900 COVID-19 samples/day
- Deliver same day, on-site results using the gold standard sensitivity and specificity of PCR technology.
- Near patients testing solutions, anywhere at anytime.

### Standard Fittings

- Mobile trailer 12ft x 6ft x 6ft 6in
- Power supply to connect to buildings including back up battery power supply for freezer
- Wifi, tracking device and alarm system
- Climate control unit
- Air sterilisation system
- -20°C laboratory freezer
- 6x q16 qPCR instruments
- Storage available for 2,600 tests and associated equipment
- Sample deposit hatch
- Externally accessible clinical waste bin compartments
- All required first aid & decontamination equipment

### Additional Extras

- Additional qPCR equipment – up to 9x q16 or 6x q32
- CL2 Biosafety hood
- Larger freezer
- Double room awning entry
- External handwash unit (8.5L capacity, hot and cold taps)
- Generator for remote locations
- Personalised branding

Contact us today for more tailored fittings to meet your testing needs!

**Call +44 (0) 02380 748 830**  
or email to  
**[enquiries@primerdesign.co.uk](mailto:enquiries@primerdesign.co.uk)**