

NOVACYT

The Choice is **YOURS!**

Introducing the PathFlow® COVID-19 Rapid Antigen Self-Test, a single-use test kit to detect SARS-CoV-2 with oral fluid sample.

So that you can choose not to swab anymore.




Rapid results in **15 minutes**



High diagnostic sensitivity of **93.0%** and diagnostic specificity of **100%** with total accuracy of **97.1%**




Easy-to-use and non-invasive



NOVACYT

COVID-19 results in 4 simple steps:



Step 1: Spit into the funnel after clearing your throat by coughing deeply

Step 2: Add the buffer into the tube

Step 3: Add 2 drops of the solution to the well labelled 'S'

Step 4: Read results at 15 minutes

Novacyt COVID-19 total solution



Novacyt can offer a total COVID-19 solution, from Instruments to self-tests



CONTACT US TODAY at sales@novacyt.com to discuss your testing needs.